AMMONIA-NITROGEN KIT

COLOR CHART, 0.1-4.0 ppm

CODE 5864-CN-01

QUANTITY	CONTENTS	CODE
50	*Ammonia #1 Tablets	*3968A-CN-H
50	*Ammonia #2 Tablets	*3969A-CN-H
2	Test Tubes, plastic, w/cap	0106
1	Ammonia-Nitrogen Color Chart, 0.1 - 4.0 ppm	6665-02-CC

***WARNING:** Reagents marked with an * are considered to be potential health hazards. To view or print a Safety Data Sheet (SDS) for these reagents go to www.lamotte.com. Search for the four digit reagent code number listed on the reagent label, in the contents list or in the test procedures. Omit any letter that follows or precedes the four digit code number. For example, if the code is 4450WT-H, search 4450. To obtain a printed copy, contact LaMotte by email, phone or fax.

Emergency information for all LaMotte reagents is available from Chem-Tel: (US, 1-800-255-3924) (International, call collect, 813-248-0585).

To order individual reagents or test kit components, use the specified code number.

PROCEDURE

1. Fill test tube (0106) to 5 mL line with sample water.
2. Add one *Ammonia #1 Tablet (3968A-CN) and one *Ammonia #2 Tablet (3969A-CN). Cap and mix until tablets disintegrate. Wait 5 minutes.
3. Hold test tube flat against the white section of the Ammonia-Nitrogen Color Chart (6665-02-CC). Match sample color to a color standard. Record as ppm Ammonia-Nitrogen.

NOTE: Sample may be turbid. This will not affect the test results. To convert result to Ammonia, multiply reading by 1.3. Record as ppm Ammonia.

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LaMOTTE COMPANY

Helping People Solve Analytical Challenges

PO Box 329 · Chestertown · Maryland · 21620 · USA
800-344-3100 · 410-778-3100 [Outside U.S.A.] · Fax 410-778-6394
www.lamotte.com
65864-CN-01 09.18