

# Quick Start | UV/VIS Spectrophotometer

1. Press the power switch on the back of the instrument to turn the instrument on. The <b>Initializing</b> screen will appear.	Initializing	15 : 00
	Booting System : Locating filter ... Warmup 15 min...	
	LAMOTTE SMART SPECTRO	

2. Press <b>ENTER</b> to select Programmed Tests.	Programmed	
	1 Sequence 1	
	2 Sequence 2	
	3 Sequence 3	
	4 All Tests	

3. Scroll to and press <b>ENTER</b> to select All Tests.	All Tests	
	1 Alkalinity-UDV	
	2 Aluminum	
	3 Ammonia-N LF	
	Press "Enter" to Run	

4. Scroll to the desired test. The spectrophotometer is ready to scan the blank. The proper wavelength has been selected.	All Tests	
	13 Ca & Mg Hard-UDV	
	14 Carbohydrazide	
	15 Chlorine	
	Press "Enter" to Run	

5. Insert the blank. Press <b>ENTER</b> to scan the blank. Wait for the instrument to blank. The blank has been stored.	Chlorine		515nm
	0.000A		99.9%T
	No.	Abs	ppm

6. Insert the reacted sample. Press <b>ENTER</b> to scan the sample. The result will be displayed.	Chlorine		515nm
	0.209		
	No.	Abs	ppm
	*01	0.212	0.309

7. Press <b>PRINT</b> to print. Insert another sample and press <b>ENTER</b> . Press <b>ESCAPE</b> to exit. Or turn spectrophotometer <b>OFF</b> .	
--	--