



Water Quality Assessment

Overview Topics

Watersheds

Everyone lives in a watershed. A watershed is an area of land where all of the water that is under the land or drains off the land flows downhill to the same place – a common waterway. Water quality is affected by how the surrounding land is used and, because water moves throughout a watershed, natural and human activities that affect the water quality upstream in the watershed can affect the water quality downstream.






Topographic Map

- Large scale map
- Elevation is shown in three dimensions
- Shows natural and manmade features
- Shows small streams as well as large water bodies
- Used to estimate watershed boundaries

You can identify your local watershed on internet websites.






Pollution

The Environmental Protection Agency defines pollution as...

“Any substance that exists in the environment that is undesirable or harmful to that environment”



The Water Quality Index (WQI)

- Developed by the EPA
- Dimensionless number
- Determines overall water quality
- Compare water quality at multiple sites
- Detect water quality trends
- Identify areas of water quality problems




Calculating the WQI

1. Nine physical, chemical, and microbiological parameters are measured:

- Dissolved Oxygen
- Fecal Coliform
- Biochemical Oxygen Demand
- pH
- Temperature Change
- Nitrate
- Phosphate
- Turbidity
- Total Dissolved Solids



2. Each test result is transferred to a weighting curve where a numerical value, or Q-value, is obtained.



A weighting factor is a mathematical value that sets the relative importance of the parameter to the overall water quality.



3. The nine resulting Q-values are used to calculate the Water Quality Index for the sampling site.

The Water Quality Index scale ranges from 1 to 100, with higher numbers indicating better water quality.

Water Quality Index Worksheet

Parameter	Test Result	Units	Q-value	Weighting Factor	Weighted Q-Value
DO	100	% Saturation	99	0.17	16.8
Fecal Coliform	500,000	CFU/100 mL	2	0.16	0.3
BOD	16.8	mg/L	15	0.11	1.7
pH	6.4	pH units	64	0.11	7.0
Temp Change	0.5	°C	90	0.10	9.0
Nitrate	9.4	mg/L NO ₃	55	0.10	5.5
Phosphate	0.4	mg/L PO ₄	72	0.10	7.2
Turbidity	24	NTU	57	0.08	4.6
TDS	218	mg/L	71	0.07	5.0

Water Quality Index = 56.9
 Water Quality Rating = Average

WQI Rating	
WQI	Rating
0-25	Poor
26-50	Fair
51-70	Average
71-90	Good
91-100	Excellent

Water Quality Assessment Curriculum Package



- Thermometer
- Dissolved Oxygen Kit
- pH Kit
- Nitrate-Nitrogen Kit

- Phosphate Kit
- Turbidity Kit
- Alkalinity Kit
- Salt/TDS/Temp Tracer

- Nutrient – TTC/MacConkey BioPaddles®
- Water Quality Index CD



Unique Feature!



Nutrient TTC/MacConkey
BioPaddles® for coliform
testing

- No Petri dishes
- No inoculating loops
- No Bunsen burners
- No refrigeration
- Faster
- SAFER



The exploration of the Water Quality Index teaches students STEM-based skills:

- Analytical skills - research a topic, develop a plan, draw conclusions
- Organization skills – record data
- Science skills - recognize cause and effect, use facts to prove a hypothesis
- Leadership skills – work as a team
- Mathematic skills – measurements and calculations
- Attention to detail - follow procedures, record data
- Technical skills - troubleshoot problems, build equipment, utilize internet sites
- Creative abilities - solve problems, design equipment
- Communication skills - interact within teams, pool class data, prepare reports


Skills will be applied to four classroom activities and one field activity...



Classroom activities use EPA data from Kansas River Sites

Activity 1	Determine the Water Quality Index for the Kansas River Watershed
Activity 2	Do Water Quality Values Change over a 24 Hour Period?
Activity 3	Do Water Quality Values Change Within a Watercourse Cross-Section?
Activity 4	Environmental Impacts on Mill Creek Tributary

1. Students review background materials





2. Students create spreadsheets with the reported data to calculate the Water Quality Index ...

KANSAS CITY - STATION 1 - 24 HOUR DATA
SS2

Date	Station	Mile Marker	Location	Time	Weather Observations	Flow	DO	Fecal Coli	BOD ₅	pH	Temp (Loc 1)	Temp (Loc 2)	Nitrate	Phosphate	Turbidity	TDS
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	2:00	Clear skies; -2°C	1980	10.8	2,000	8.3	7.7	9.5	9.5	2.000	0.56	20	450
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	4:00	Clear skies; -2°C	1980	10.8	2,000	7.6	7.7	9.5	9.5	2.000	0.43	17	435
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	6:00	Clear skies; 0°C	1980	10.6	1,000	8.9	7.8	9.5	9.5	2.000	0.80	15	460
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	8:00	Clear skies; 1°C	1970	11.0	3,000	8.9	8.0	9.0	9.0	2.000	0.43	13	475
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	10:00	Light Rain; 3°C	1940	10.1	5,000	8.3	7.8	9.0	9.0	2.000	0.61	22	435
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	12:00	Light Rain; 4°C	1940	11.1	2,500	7.6	8.0	9.0	9.0	2.000	0.13	51	425
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	14:00	Drizzle; 5°C	1960	11.1	8,000	6.4	7.7	9.0	9.0	2.000	0.62	38	458
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	16:00	Clear skies; 6°C	1980	11.1	5,000	6.4	8.0	9.0	9.0	2.000	0.46	30	440
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	18:00	Clear skies; 5°C	1930	10.7	2,500	7.3	8.0	8.5	8.5	2.000	0.59	34	435
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	20:00	Clear skies; 3°C	1930	13.8	5,000	8.0	7.8	9.0	9.0	2.000	0.44	25	470
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	22:00	Clear skies; 2°C	1900	10.4	1,000	7.9	8.0	9.0	9.0	2.000	0.55	25	425
6-Nov	1	0	Kansas City Latitude 39°05' Longitude 94°36'	0:00	Clear skies	1910	10.2	5,000	8.3	7.8	9.0	9.0	2.000	0.43	25	440




...and analyze the data and create graphs to answer questions such as:

- Based on the sampling data, can any conclusions be drawn regarding sampling time?
- Does the data show a significant difference in water quality over time?
- Based on the reported data, can rainfall influence the Water Quality Index?
- What water quality target parameters fell outside of the target values?

And discuss:

- What are the possible reasons for water quality variances?
- What would be a possible plan for future sampling that might yield more significant data?
- How would you design a study that would provide additional guidance regarding river and stream sampling points?

Completed spreadsheets and answers to questions are provided in the Teacher Guide



TEACHERS NOTE:
You may wish to assign data from specific dates to student groups and have groups share the Water Quality Index values that they calculate to create a composite spreadsheet.

4. Create a Statistical Summary Spreadsheet
Have students create a spreadsheet that presents a statistical summary of the watershed WQI values, including:

Sample size (n)	n = 10
Mean (\bar{x})	35.16
Standard deviation (SD)	3.08
Standard error of the mean (SEM)	1.0
Data range (max / min)	39.04 - 29.84

5. Summarize the WQI
Have students summarize the WQI for each sampling date.

Sampling Date	WQI
9-26-11	36.08
10-4-11	37.39
10-12-11	37.55
10-20-11	34.10
11-4-11	38.23
11-5-11	34.78
11-13-11	39.05
11-21-11	33.56
11-29-11	31.00
12-7-11	29.84

Have students calculate the ave of September thru early Decem
WQI = 35.16

Have students determine whethr

Parameter	Tar
Dissolved Oxygen	9.5
pH	6.5
BOD	< 1
Nitrate	4 m
Phosphate	0.1

Have students discuss whether source.
There is a feedlot located appr is occurring.

KANSAS CITY - STATION 1 - 24 HOUR DATA
SS2 - ANSWER KEY

Date	Station	Location	Time	Weather	Flow (cfs)	DO (mg/L)	DO % Saturation	Q	Fecal Coll (CFU/100 ml)	Q	BOD ₅ (mg/L)	Q	pH	Q	Temp (°C) (Surf)	Temp (°C) (Surf)	Temp (°C) (Surf)	Q	Nitrate (mg/L NO ₃ -N)	Q	Phosphate (mg/L PO ₄ -P)	Q	Turbidity (FTU)	Q	TDS (ppm)	Q	WQI	WQI Rating
6-Nov	1	Kansas City Lumberton St SW	2:00	Clear skies, 2°C	1980	19.8	88	92.40	19.00	8.3	41.14	7.7	89.80	9.5	9.5	0.0	92.00	2.000	88.40	0.56	85.96	20	60.00	450	39.00	85.16	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	4:00	Clear skies, 2°C	1980	19.8	85	90.75	18.00	7.6	44.08	7.7	89.80	9.5	9.5	0.0	92.00	2.000	88.40	0.43	73.63	17	65.40	435	41.40	86.57	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	6:00	Clear skies, 3°C	1980	19.8	85	90.75	22.00	8.8	38.62	7.8	88.20	9.5	9.5	0.0	92.00	2.000	88.40	0.80	51.80	15	69.00	460	37.40	84.26	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	8:00	Clear skies, 1°C	1970	19.7	82	94.60	15.79	8.8	38.62	8.0	88.00	9.0	9.0	0.0	92.00	2.000	88.40	0.43	73.63	13	72.60	475	35.00	83.87	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	10:00	Light Rain, 3°C	1940	19.4	73	77.50	13.00	8.3	41.14	7.8	88.20	9.0	9.0	0.0	92.00	2.000	88.40	0.61	63.01	22	58.60	435	41.40	81.41	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	12:00	Light Rain, 4°C	1940	19.4	82	85.10	16.78	7.8	44.08	8.0	88.00	9.0	9.0	0.0	92.00	2.000	88.40	0.13	91.33	37	60.00	425	43.00	85.22	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	14:00	Drizzle, 5°C	1960	19.6	73	77.50	10.97	6.4	43.12	7.7	89.80	9.0	9.0	0.0	92.00	2.000	88.40	0.62	62.42	38	47.40	458	37.72	80.93	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	16:00	Clear skies, 3°C	1980	19.8	95	96.25	13.00	6.4	43.12	8.0	88.00	9.0	9.0	0.0	92.00	2.000	88.40	0.48	71.86	30	63.00	440	40.00	85.50	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	18:00	Clear skies, 2°C	1930	19.3	95	96.80	16.78	7.3	45.34	8.0	88.00	9.5	9.5	0.0	92.00	2.000	88.40	0.59	64.19	34	50.20	435	41.40	84.83	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	20:00	Clear skies, 2°C	1930	19.3	105	96.75	13.00	8.0	42.40	7.8	88.20	9.0	9.0	0.0	92.00	2.000	88.40	0.44	75.04	28	56.50	470	35.80	83.28	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	22:00	Clear skies, 2°C	1900	19.0	105	96.75	1.000	7.8	42.82	8.0	88.00	9.0	9.0	0.0	92.00	2.000	88.40	0.55	66.55	25	56.50	425	43.00	88.23	AVERAGE	
6-Nov	1	Kansas City Lumberton St SW	0:00	Clear skies	1910	19.2	105	96.75	5.000	13.00	8.3	41.14	7.8	88.20	9.0	9.0	0.0	92.00	2.000	88.40	0.43	75.63	25	56.50	455	40.00	85.52	AVERAGE
Mean	1950	11.0	90	91.33	2.000	16.80	7.8	43.14	7.8	87.77	8.1	8.1	8.0	9.0	9.0	0.0	92.00	2.000	88.40	0.50	68.20	26	56.50	446	38.89	84.73	AVERAGE	
SD	29	1.0	11	7.02	0.000	0.00	0.8	3.46	0.11	2.10	0.3	0.3	0.0	0.00	0.000	0.00	0.16	8.55	11	9.54	17	5.64	17	5.64	1.75			
Min	6	0	5	0.0	400	1.0	0	1.0	0	0.0	0	0	0	0.0	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0	0.00	0	0.00	0.00		
Max	1960	10.1	75	77.50	1.000	18.87	8.4	38.62	7.7	44.08	8.4	8.5	8.0	9.0	9.0	0.0	92.00	2.000	88.40	0.80	51.80	15	69.00	425	35.00	83.88		
Max	1980	13.8	105	96.80	0.000	22.00	8.8	49.12	8.0	89.80	9.5	9.5	9.0	9.0	9.0	0.0	92.00	2.000	88.40	0.80	51.80	15	72.60	475	43.00	86.57		



Activity 5


Investigate Water Quality in a Local Watershed

Field Activity



Identify a Local Watershed


Students do a search or use links to internet resources to identify a local watershed and determine the local watershed boundary.






Student explore internet sites to research more about the site:

- Land Use (GIS)
- Water quality data from citizen monitoring groups
- Weather events
- Water quality data for nearby monitoring stations





Students choose a sampling site,





... list possible point and nonpoint pollution sources,...

- Urban runoff
- Sewage
- Industrial wastewater
- Thermal pollution
- Debris from upstream
- Secondary-treated effluent
- Waste discharges



...and use satellite photos from internet sites to document the sampling location.



Collect a Water Sample

Students collect a sample at the site



...or a sample can be brought to the classroom.





**Change in
Temperature**
0.5°C
Q=90



**Biochemical
Oxygen Demand**
16.8 mg/L
Q=15



Coliform
500,000 CFU/100 mL
Q=2



Students enter test results on the Field Data Sheet and the Water Quality Index Worksheet.

Water Quality Index Worksheet

Parameter	Test Result	Units	Q-value	Weighting Factor	Weighted Q-Value
DO	100	% Saturation	99	0.17	16.8
Fecal Coliform	500,000	CFU/100 mL	2	0.16	0.3
BOD	16.8	mg/L	15	0.11	1.7
pH		pH units		0.11	
Temp Change	0.5	°C	90	0.10	9.0
Nitrate		mg/L NO ₃	53	0.10	
Phosphate		mg/L PO ₄	72	0.10	
Turbidity		NTU	57	0.08	
TDS		mg/L	71	0.07	

Water Quality Index = _____

Water Quality Rating = _____

WQI Rating	
WQI	Rating
0-25	Poor
26-50	Fair
51-70	Average
71-90	Good
91-100	Excellent



MISSING DATA!

BECOME THE STUDENTS!

TEST THE WATER SAMPLES FOR THE REMAINING TEST FACTORS.

Water Quality Index Worksheet

Parameter	Test Result	Units	Q-value	Weighting Factor	Weighted Q-Value
DO	100	% Saturation	75	0.17	16.8
Fecal Coliform	500,000	CFU/100 mL	2	0.16	0.3
BOD	16.8	mg/L	15	0.11	1.7
pH		pH units		0.11	
Temp Change	0.5	°C	90	0.10	9.0
Nitrate		mg/L NO ₃		0.10	
Phosphate		mg/L PO ₄		0.10	
Turbidity		NTU		0.08	
TDS		mg/L		0.07	

Water Quality Index =

Water Quality Rating =

WQI Rating	
WQI	Rating
0-25	Poor
26-50	Fair
51-70	Average
71-90	Good
91-100	Excellent



1. Follow the Field Test Procedure to test the water sample from the site.
2. Enter the test result on the Field Data Sheet.

MEASURE NITRATE

For Use with: LaMotte Nitrate Nitrogen Tablet Kit (Code 3354-01). See Resources for links to complete test kit instructions and chemical reactions

Nitrogen is required by all living animals and plants to build protein. Aqi or they eat other aquatic organisms that consume aquatic plants. Nitro decomposition of plants and animals or the excrement of living animals amount of nitrogen in areas where they are plentiful.

Nitrogen is a nutrient that acts as a fertilizer for aquatic plants and, along causing water quality problems.

Drinking water containing high nitrate levels can affect the ability of blood animals are especially susceptible.

Nitrogen enters the water from improperly functioning septic systems, a of fertilizer from lawns and crops.

Unpolluted waters generally have a nitrate-nitrogen level below 1 ppm. I are considered unsafe for drinking water.

Nitrogen occurs in natural waters as nitrate (NO_3^-), nitrite (NO_2^-). Ammonia results are expressed as nitrate-nitrogen ($\text{NO}_3\text{-N}$), meaning "nitrogen the expressions refer to the same chemical but express the concentration in

MATERIALS

Nitrate Nitrogen Tablet Kit	Code 3354-01
Water collection bottle, at least 1 liter	
Thermometer	Code 1066
Field Data Form	
Water Sampler (optional)	Code 1054-DO

SAFE

-
-
-
-

WARNING! This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully. Not to be used by children except under adult supervision.

MEASURE TURBIDITY (TOTAL SUSPENDED SOLIDS-TSS)

For Use with: LaMotte Turbidity Kit (Code 7519-01) or LaMotte Secchi Disk (Code 0171-CL)(alternative method). See Resources for link to complete test kit instructions.

Turbidity is a measurement of the relative clarity of water. Turbidity is caused when suspended solids scatter light as it passes through water. Turbidity should not be confused with color, since darkly colored water can still be clear and not turbid.

In addition to blocking out the light needed by submerged aquatic vegetation, burying eggs and bottom dwelling creatures, suspended sediment can carry nutrients and pesticides throughout the water system. High levels of turbidity can damage gills and interfere with the ability of fish to see food. Suspended particles near the water surface absorb additional heat from sunlight, raising surface water temperature.

There are many possible sources of turbidity. Sediment from disturbed or corroded soil can cloud water and microscopic plankton contribute to high turbidity levels when their numbers increase due to excess nutrients and sunlight. Waste discharge, urban runoff, and large quantities of bottom feeders or boat traffic that stir up bottom sediment can increase turbidity levels.

Turbidity can be measured in several ways. One method uses a turbidimeter-- an instrument that electronically measures how light is scattered when it passes through a water sample. The more light that is scattered, the greater the turbidity reading will be.

WARNING! This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully. Not to be used by children except under adult supervision.

The Secchi Disk provides a very simple means of making transparency determinations in natural waters. The Secchi Disk is a black and white circular plastic plate, 20 cm in diameter that is divided into four pie-shaped quadrants that are alternately colored black and white. The Secchi disk is lowered into the water until it can no longer be seen. The depth at which it vanishes is measured and recorded. Secchi disks do not work well in shallow waters or in waters where the current is too strong.

For shallower waters, a disappearing dot method can be used. The water sample is added to a tube with a target at the bottom. A calibrated turbidity standard is added to clear water in a second tube until the observed "fuzziness" of both targets is the same.

Moderately low levels of turbidity may indicate a healthy, well-functioning ecosystem in which plankton flourish at a reasonable level to form the foundation of the food web. High turbidity is an indicator of either runoff from disturbed or eroded soil or blooms of microscopic organisms due to high nutrient inputs. Very clear water is typical of the open ocean and oligotrophic lakes, supporting only sparse plant and animal life.

Turbidity levels in drinking water should be less than 0.5 JTU and typical groundwater levels are less than 1.0 JTU.

FIELD COLLECTION

(If using the Secchi disk, go to Alternative Field Collection)

MATERIALS

Turbidity Test Kit	Code 7519-01
Water collection bottle, at least 1 liter	
Thermometer	Code 1066
Field Data Form	
Water Sampler (optional)	Code 1054-DO

SAFETY

- Wear splash-certified safety goggles and gloves.
- Read the labels on all reagents
- Wash hands before and after contacting water samples and performing tests. Poor hygiene is a common reason samples become contaminated.
- In the event of an accident or suspected poisoning, immediately call the Poison Center phone number in the front of your local telephone directory or call a physician. Additional information for all LaMotte reagents is available in the United States, Canada, Puerto Rico, and the US Virgin Islands from Chem-Tel by calling 1-800-255-3924. For other areas, call 813-248-0585 collect to contact Chem-Tel's International access number. Each reagent can be identified by the four digit number listed on the upper left corner of the reagent label, in the contents list and in the test procedures.



Calculate the Water Quality Index

1. Pool group data to complete the Water Quality Index Worksheet.
2. Set up a spreadsheet, use the Water Quality Index Calculator, or use the Water Quality Index Worksheet calculate the Water Quality Index.

Water Quality Index Calculator

Parameter	TEST RESULT	Units	Q-value	Weight	Weighted Q-value
DO		% Saturation	NR	0.17	NR
Fecal Coliform		CFU/100 mL	NR	0.16	NR
BOD		mg/L	NR	0.11	NR
pH		pH units	NR	0.11	NR
Temp Change		°C	NR	0.10	NR
Nitrogen		mg/L NO ₃	NR	0.10	NR
Phosphate		mg/L PO ₄	NR	0.10	NR
Turbidity		JTU or in/ft	NR	0.08	NR
TDS		ppm	NR	0.07	NR

NR = Not Reported

WEIGHTS: 0 0.00

Water Quality Index = NR

Water Quality Rating = NR



WHAT IS THE WATER QUALITY OF THE LOS ANGELES RIVER BETWEEN WARDLOW ROAD AND WILLOW ROAD?

WHAT COULD HAVE INFLUENCED THE WATER QUALITY?




Rainfall?

Restricted water flow?

Sewer and industrial
wastewater discharge?

Urban runoff?

When you drive over the Los Angeles River and through the city of Long Beach will you look at the area the same way as you would have 90 minutes ago?





The concentrations of water sample parameters, physical data, and pollution sources used in the example activity were based on actual data that was collected from the Los Angeles River between Wardlow Road and Willow Street on October 11, 2012





AP® The Water Quality
Assessment Curriculum Module
Order Code 5845



AP® The Water Quality
Assessment Curriculum Package
Order Code 5845-PKG



www.lamotte.com

Booth #716