



QUANTITY	CONTENTS	CODE	
15 mL	Iodine Reagent #1	7254-E	
15 mL	*Iodine Reagent #2	*7255-E	
60 mL	Iodine Reagent #3	6406PS-H	
1	Test Tube, 5-10-15-20-25-30 mL, plastic, w/cap	0715	<p>*Reagent is a potential health hazard. READ SDS: lamotte.com</p> <p>Emergency information:</p> <p>Chem-Tel USA 1-800-255-3924</p> <p>Int'l, call collect, 813-248-0585</p> <div><div>READ SDS</div><div></div><div></div></div>
1	Pipet, dropping, plastic	0352	
To order individual reagents or test kit components, use the specified code number.			

Warning! This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully.
Not to be used by children except under adult supervision.

PROCEDURE

1. Fill test tube [0715] to 25 mL line with sample water.
 2. Add 10 drops of Iodine Reagent #1 [7254] and 10 drops of *Iodine Reagent #2 [7255]. Cap and swirl to mix. Solution will turn blue.
 3. Fill pipet [0352] with Iodine Reagent #3 [6406PS]. While gently swirling test tube, add Iodine Reagent #3 one drop at a time, until blue color disappears. Count the number of drops added. Hold pipet vertically.
 4. Multiply number of drops used in Step 3 by 2.5. Record as ppm Iodine.
- Number of Drops x 2.5 = ppm Iodine**