

TSA/RB

Code 5552

Tryptic Soy Agar (TSA)

Rose Bengal Chloramphenicol (RB)

USE:

Cultivation of a wide variety of aerobic and anaerobic microorganisms¹. (TSA)

Selective enumeration and cultivation of yeasts, molds and Actinomycetes from food and other surfaces. (RB)

Side 1: Tryptic Soy Agar (TSA) (off white, opaque)

(Side 1 of each paddle is marked with an indented laser line)

Side 2: Rose Bengal Chloramphenicol Agar (RB) (pink)



APPLICATION

Tryptic Soy Agar is recommended in multiple water & wastewater applications², and numerous standard methods for food testing.³ TSA is commonly used as a maintenance medium for culture collections and testing bacterial contaminants in cosmetics.⁴ Rose Bengal Chloramphenicol Agar is recommended in Standard Methods for the Examination of Water and Wastewater for the enumeration of yeasts and molds from foods and water.

PADDLE AGARS

Tryptic Soy Agar (TSA) – In 1955, Leavitt et al.⁵ discovered that Tryptic Soy Agar facilitated vigorous growth of aerobic and anaerobic microorganisms. This is an enriched media, suitable to support fastidious heterotrophs.

Rose Bengal Chloramphenicol Agar (RB) – Selective medium for the enumeration of fungi. This formula is prepared with a neutral pH and supplemented with chloramphenicol as the selective agent in fungal medium. RB Agar is also referred to as Rose Bengal Agar and Rose Bengal-Malt Extract Agar. Agar and a proprietary polymer are the solidifying agents.

¹ United States Pharmacopeial Convention. 1995. The United States Pharmacopeia, 23rd ed. The United State Pharmacopeial Convention, Rockville, MD.

² Greenberg, A. E., L. S. Clesceri, and A. D. Eaton (eds.). 1995. Standard methods for the examination of water and wastewater, 19th ed. American Public Health Association, Washington, D.C.

³ U.S. Food and Drug Administration. Bacteriological analytical manual, 8th ed., AOAC International, Gaithersburg, MD.

⁴ Curry, A. S., G. G. Joyce, and G. N. McEwen, Jr. 1993. CTFA Microbiology guidelines. The Cosmetic, Toiletry, and Fragrance Association, Inc. Washington, D.C.

⁵ Leavitt, J. M., I. J. Naidorf and P. Shugaevsky. 1955. The undetected anaerobe in endodontics: a sensitive medium for detection of both aerobes and anaerobes. The NY J. Dentist. 25:377-382.

CULTURE CONTROLS

10-100 inoculum (CFU)

	TSA Agar	RB Agar
<i>Aspergillus niger</i>	GROWTH	GROWTH
<i>Bacillus subtilis</i>	GROWTH	INHIBITED
<i>Candida albicans</i>	GROWTH	INHIBITED
<i>Escherichia coli</i>	GROWTH	INHIBITED
<i>Pseudomonas aeruginosa</i>	GROWTH	INHIBITED
<i>Staphylococcus aureus</i>	GROWTH	INHIBITED

STORAGE / EXPIRATION

Store tightly sealed BioPaddles® in a cool, dry location. Shield from direct sunlight. Store BioPaddles® at room temperature (65 - 77°F/18 - 25°C). Avoid sudden temperature changes. Temperature fluctuations may result in condensation settling at the bottom of the vial. This will not affect the culture properties but could reduce the shelf-life or cause the agar to separate from the plastic paddle support. Do not refrigerate or store at temperatures above 80°F/27°C. Refrigeration may result in water condensation. Avoid freezing. Freezing can promote excess water loss and variation in media surface due to crystal formation. If freezing occurs, wrap BioPaddle in vial in thick towel and thaw at room temperature for 3-6 hours.

Refer to Best Before End date (See: BBE stamped on vial). Discard if paddle agar appears oxidized and darker than the expected color or if contaminants appear. The expiration date is based on medium in an intact container that is stored as directed.

SAMPLING

Liquids: Twist to remove paddle from vial. Fill vial to 40 mL fill line with the liquid to be sampled. The 40 mL volume can be used to calculate Total Viable Count (TVC) and/or Total Colony Count (TCC). Replace paddle. Allow a contact time of 15 seconds. Remove the paddle. Empty the vial. Replace the paddle in the vial.



Surfaces: Recovery rate is about 50%. Twist paddle to remove from vial. To ensure an accurate recovery, touch the paddle surface (10 cm²) to the test surface twice to cover a 20 cm² area (2 X 10 cm²). Allow 15 second contact time. Replace paddle in vial.


INCUBATION

Temperature	Minimum Period	Optimal Period
35°C	18 hours	24 hours
20-25°C	5 days	7 days






Incubate at 30 - 35°C for 18 hours. Growth will be present on the TSA agar, as it is enriched and non-selective. Growth may or may not be present on the Rose Bengal agar after this brief period as it is selective for yeasts and molds which tend to like cooler temperatures and grow more slowly in general. For this reason, re-incubate and check again after 2 to 4 days.

COLONY MEASURING

Each BioPaddles® paddle has molded media attachment points that are 4mm in length (point-to-point). This feature provides a useful guidepost to estimating nearby colony size.















IDENTIFICATION

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Actinomyces bovis</i>	<ul style="list-style-type: none"> • Lactose (+ / v) • Indole (v) ◆ • Oxidase (-) ◆ • Catalase (-) ◆ • Urease (-) ◆ • Gram (+) Rod 	PARTIAL - COMPLETE INHIBITION	---	---	++	<ul style="list-style-type: none"> • Opaque/tan-gray • Convex • Glossy • Entire • 1-3mm 	
<i>Alternaria spp.</i>	<ul style="list-style-type: none"> • Catalase (+) ◆ • Ascomycete 	+	<ul style="list-style-type: none"> • Downy to wooly; flat • Grayish, short, aerial hyphae • later becomes greenish black or olive brown with a light border • 3-9cm 		++	<ul style="list-style-type: none"> • Downy to wooly; flat • Grayish, short, aerial hyphae • later becomes greenish black or olive brown with a light border • 3-9cm 	
<i>Aspergillus niger</i>	<ul style="list-style-type: none"> • Catalase (+) ◆ • Ascomycete 	+	<ul style="list-style-type: none"> • Granular • Jet black conidia w/yellow/gray hyphae • 3-9cm 		+++	<ul style="list-style-type: none"> • Granular • White, w jet black fruiting bodies w/ yellow/gray hyphae • 3-9cm 	

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles®. All rights reserved.










LaMotte_BioPaddles_TSA_RB 07.22

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Aspergillus flavus</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Granular to wooly Yellow, yellow-green or yellow-brown pigment 3-9cm 		+++	<ul style="list-style-type: none"> Granular to wooly Yellow, yellow-green or yellow-brown pigment 3-9cm 	
<i>Aspergillus fumigatus</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Granular to cottony Blue-green, green-gray, green-brown pigment 3-9cm 		+++	<ul style="list-style-type: none"> Granular to cottony Blue-green, green-gray, green-brown pigment 3-9cm 	
<i>Aspergillus terreus</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Granular, radially rugose (wrinkled) Cinnamon buff, brown pigment 3-9cm 		+++	<ul style="list-style-type: none"> Granular, radially rugose (wrinkled) Cinnamon buff, brown pigment 3-9cm 	
<i>Botrytis spp.</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Wooly white, grey/brown pigment 3-9cm 		+++	<ul style="list-style-type: none"> Wooly white, grey/brown pigment 3-9cm 	
<i>Candida albicans</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+++	<ul style="list-style-type: none"> Cream Convex Glossy Entire 1-2mm 		+++	<ul style="list-style-type: none"> Pink Spreading 6mm 	
<i>Chaetomium spp.</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Wooly white, grey/olive pigment 3-5cm 		+++	<ul style="list-style-type: none"> Wooly white, grey/olive pigment 3-5cm 	

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles®. All rights reserved.













LaMotte_BioPaddles_TSA_RB 07.22

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Cladosporium</i> spp.	<ul style="list-style-type: none"> • Catalase (+) ◆ • Ascomycete 	+	<ul style="list-style-type: none"> • Slow-growing • Mostly olivaceous-brown to blackish but also sometimes grey, buff or brown, suede-like to floccise (wooly), often becoming powdery due to the production of abundant conidia • Vegetative hyphae, coniciophores and conidia are equally pigmented 		+	<ul style="list-style-type: none"> • Slow-growing • Mostly olivaceous-brown to blackish but also sometimes grey, buff or brown, suede-like to floccise (wooly), often becoming powdery due to the production of abundant conidia • Vegetative hyphae, coniciophores and conidia are equally pigmented 	
<i>Epicoccum</i> spp.	<ul style="list-style-type: none"> • Catalase (+) ◆ • Ascomycete 	+++	<ul style="list-style-type: none"> • Woolly/Cottony/Felty • yellow, orange, red, brown pigment • 3-5+cm 		+++	<ul style="list-style-type: none"> • Woolly/Cottony/Felty • Purplish / black • 3-5+cm 	
<i>Escherichia coli</i>	<ul style="list-style-type: none"> • Lactose (+) • Indole (+) ◆ • Oxidase (-) ◆ • Catalase (+) ◆ • Urease (-) ◆ • Gram (-) Rod 	+++	<ul style="list-style-type: none"> • Translucent/translucent; spreading • Convex, glossy, moist • Entire • 1-2 mm 		INHIBITED	---	---
<i>Enterobacter aerogenes</i>	<ul style="list-style-type: none"> • Lactose (+) • Indole (-) ◆ • Oxidase (-) ◆ • Catalase (+) ◆ • Urease (-) ◆ • Gram (-) Rod 	+++	<ul style="list-style-type: none"> • Clear to amber • Convex • Glossy • Entire • 2-4mm 		++	<ul style="list-style-type: none"> • Pink to Red • Convex • Glossy • Entire • 2-4mm 	
<i>Fusarium</i> spp.	<ul style="list-style-type: none"> • Catalase (+) ◆ • Ascomycete 	+++	<ul style="list-style-type: none"> • Fast-growing (4.5 cm in 4 days) • Wooly, white to yellow, pink, red or purple shades • 4+ cm 		+++	<ul style="list-style-type: none"> • Fast-growing (4.5 cm in 4 days) • Wooly, white to yellow, pink, red or purple shades • 4+ cm 	

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles®. All rights reserved.










LaMotte_BioPaddles_TSA_RB 07.22

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Microsporium</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Glabrous (smooth)/Downy/Wooly/Powdery White at first, later becoming grayish yellow to blue green with age Wrinkled w/ age 1-9cm 		+	<ul style="list-style-type: none"> Glabrous (smooth)/Downy/Wooly/Powdery White at first, later becoming grayish yellow to blue green with age Wrinkled w/ age 1-9cm 	
<i>Muccor</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Zygomycete 	+	<ul style="list-style-type: none"> Wooly, fast-growing Initially white, the white-yellow to various shades of grey to green showing lollipop sporangia (60-300 µm) 3-9+ cm 		+	<ul style="list-style-type: none"> Wooly, fast-growing Initially white, the white-yellow to various shades of grey to green showing lollipop sporangia (60-300 µm) 3-9+ cm 	
<i>Penicillium chrysogenum</i> (notatum)	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Granular/wooly, flat Initially white, then various shades of green, blue-green, or white-green pigment 3-9+ cm 		++	<ul style="list-style-type: none"> Granular/wooly, flat Initially white, then various shades of green, blue-green, or white-green pigment 3-9+ cm 	
<i>Penicillium roqueforti</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Granular, velvet-like, flat Initially white, then various shades of green, blue-green pigment 3-9+ cm 		++	<ul style="list-style-type: none"> Granular, velvet-like, flat Initially white, then various shades of green, blue-green pigment 3-9+ cm 	
<i>Penicillium digitatum</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+++	<ul style="list-style-type: none"> Suede-like, wooly White at first, later becoming green with age 3-9cm 		+++	<ul style="list-style-type: none"> Suede-like, wooly White at first, later becoming green with age 3-9cm 	
<i>Pithomyces</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Powdery Pale/dark grey or brown pigment 3-9cm 		+++	<ul style="list-style-type: none"> Powdery Pale/dark grey or brown pigment 3-9cm 	

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles®. All rights reserved.









LaMotte_BioPaddles_TSA_RB 07.22

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Pseudomonas aeruginosa</i>	<ul style="list-style-type: none"> Lactose (-) Indole (-) ◆ Oxidase (+) ◆ Catalase (+) ◆ Urease (-) ◆ Gram (-) Rod Fluoresces blue under long-wave UV light (400-nm) 	+++	<ul style="list-style-type: none"> Transparent to amber Circular to irregular, spreading Slightly convex, glossy Entire 0.2-0.5 mm 		INHIBITED	---	---
<i>Pseudomonas fluorescens</i>	<ul style="list-style-type: none"> Lactose (-) Indole (-) ◆ Oxidase (+) ◆ Catalase (+) ◆ Urease (-) ◆ Gram (-) Rod Fluoresces blue-green under long-wave UV light (400-nm) 	+++	<ul style="list-style-type: none"> Dull amber to yellow Irregular Umbonate, glossy, butyrous Undulate 4-6 mm 		+/-	<ul style="list-style-type: none"> Red, pink Irregular; Spreading to confluent 2-4mm 	
<i>Rhizopus spp.</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Zygomycete 	+++	<ul style="list-style-type: none"> Cottony white to blackish grey (black fruiting bodies) 3-9cm 		+++	<ul style="list-style-type: none"> Cottony white to blackish grey (black fruiting bodies) 3-9cm 	
<i>Saccharomyces cerevisiae</i>	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+++	<ul style="list-style-type: none"> Transparent to white to cream Punctiform Convex, glossy Entire 0.1-0.5 mm 		+++	<ul style="list-style-type: none"> White to pink Punctiform/circular Convex, dull Entire 0.1-0.5 mm 	
<i>Salmonella typhimurium</i>	<ul style="list-style-type: none"> Lactose (-) Indole (-) ◆ Oxidase (-) ◆ Catalase (+) ◆ Urease (-) ◆ Gram (-) Rod 	+++	<ul style="list-style-type: none"> Transparent Circular Semi-convex Entire 0.5 - 1.0mm 		INHIBITED	---	---
<i>Salmonella enteritidis</i>	<ul style="list-style-type: none"> Lactose (-) Indole (-) ◆ Oxidase (-) ◆ Catalase (+) ◆ Urease (-) ◆ Gram (-) Rod 	+	<ul style="list-style-type: none"> Transparent Circular Semi-convex Entire 0.5 - 1.0mm 		INHIBITED	---	---

For *in vitro* diagnostic use only. This product should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

©LaMotte BioPaddles®. All rights reserved.

LaMotte_BioPaddles_TSA_RB 07.22

ORGANISM		TSA			RB		
ORGANISM	PHYSIOLOGY ◆ Precision Test Strip Available	GROWTH	COLONY	IMAGE	GROWTH	COLONY	IMAGE
<i>Stachybotrys</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Wooly black (sometimes: white, pink, orange) 3-9cm 		++	<ul style="list-style-type: none"> Wooly black (sometimes: white, pink, orange) 3-9cm 	
<i>Torula</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> White Punctiform Convex, glossy Entire 0.1-0.5 mm 		+	<ul style="list-style-type: none"> White Punctiform/circular Raised Entire to irregular 0.1-0.5 mm 	
<i>Trichoderma</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	++	<ul style="list-style-type: none"> Cottony White/late scattered green or yellow-green patches (rings) 3-9cm 		++	<ul style="list-style-type: none"> Cottony White/late scattered green or yellow-green patches (rings) 3-9cm 	
<i>Trichophyton</i> spp.	<ul style="list-style-type: none"> Catalase (+) ◆ Ascomycete 	+	<ul style="list-style-type: none"> Wooly with indented borders White to brownish-tan pigment 3-9cm 		+	<ul style="list-style-type: none"> Wooly with indented borders White to brownish-tan pigment 3-9cm 	
+++ = very rich, luxurious growth expected ++ = grows + = grows slightly +/- = may grow; may be inhibited							

DISPOSAL

Twist to remove paddle from vial. Fill vial to 40 mL fill line with 1:9 dilution of household bleach (5.25% sodium hypochlorite). Replace paddle in vial. Allow 15 minute contact time. Remove paddle. Discard bleach solution. Replace paddle in vial and dispose. Alternatively, loosen cap and microwave for 30 seconds, autoclave, or incinerate.

GLOSSARY:

Catalase Test	Catalase enzyme will react with hydrogen peroxide to produce oxygen if the bacteria is catalase positive.
Lactose Test	Lactose positive bacteria can ferment available lactose in the agar producing an acid which lowers the pH. Lactose negative bacteria are non-fermenting.
Indole Test	Biochemical test to determine the ability of an organism to split indole from the amino acid tryptophan. <i>P. vulgaris</i> is indole positive while <i>P. mirabilis</i> is indole negative.
Oxidase Test	Oxidase positive bacteria contain cytochrome c oxidase which will turn an indicator dark blue. In contact with oxidase negative bacteria, the indicator will remain colorless.
Urease Test	Bacteria containing urease will hydrolyze urea to ammonia and carbon dioxide causing an alkaline environment which changes the color of a pH indicator from yellow to fuchsia.
β-D-Glucoronidase Reaction	The presence of <i>E. coli</i> is determined when both β-D-Glucoronidase and Indole are positive, and the organism is gram negative.
Gram Staining	A method for differentiating bacteria into two groups – gram positive and gram negative – based on the chemical and physical properties of their cell walls. Often the first step in identifying bacteria.